

Downloadable checklist - From Excel to LIMS: digitizing water samples, results and reports

Audience: Laboratories, QA/QC, distributors, utilities, multi-site operations · Region: Global

Prepared for water microbiology, traceability and product selection conversations.

Scope

- Audience: Laboratories, QA/QC, distributors, utilities, multi-site operations
- Region: Global

Controls to verify

- ☐ Inventory spreadsheets, reports, templates and critical data.
Status: Responsible person: Notes:
- ☐ Map the real flow from request to report and invoice.
Status: Responsible person: Notes:
- ☐ Define mandatory fields for sample, point, method, operator and lot.
Status: Responsible person: Notes:
- ☐ Separate result entry, technical review and report release.
Status: Responsible person: Notes:
- ☐ Connect CoA, customer portal and support communications.
Status: Responsible person: Notes:
- ☐ Plan migration by phases with a minimum viable workflow.
Status: Responsible person: Notes:

Recommended links

- Request technical recommendation
- View AquaVerify products
- Request AquaVerify Cloud demo

[Request technical recommendation](#)

[View AquaVerify products](#)

[Request AquaVerify Cloud demo](#)

Suggested next step

Use the completed checklist to prepare a product, platform, distributor or OEM discussion with AquaVerify.

HTML resource URL:

https://aquaverify.com/en/contact?intent=technical_recommendation&resource=excel-to-lims

This document is technical orientation material. Method validation, accreditation, competent authority requirements and the quality system of each organization remain decisive.